

Baymedix® A Adhesive and Platilon® Film Combinations Raise Performance and Efficiency.

Baymedix® A and Platilon® polyurethanes are designed for modern wound care and NPWT applications. In combination, these materials show excellent breathability, an important feature in moisture management of advanced wound dressings. The high moisture vapor transmission rate removes the need for perforation of the adhesive layer. Direct coating on film eliminates processing or lamination steps in (conversion) production.

Baymedix®

Key Features of Adhesives made from Baymedix® A:

- Adhesive peel strength easily adjusted by component ratio
- Low trauma removal without sticking to hair

Key Features of Baymedix® A and Platilon® Film:

- Excellent breathability
- Enhanced processing efficiency
- Matt surface (PE facing) of Platilon® Film



Two-Component Adhesive

Property	Baymedix® AP501
NCO content	12.3–13.3%
Viscosity at 23°C	3,300–4,700 mPa · s
Solvent	no
Monomer content	≤ 0.39%

Property	Baymedix® AR602
OH number	32.5–37.5 mg KOH/g
Viscosity at 25°C	900–1,100 mPa · s
Solvent	no
Acid number	≤ 0.04 mg KOH/g


Polyurethane Film

Property	VPT 9101 T
MVTR	~1,900 g/(day · m ²)
Tensile stress at break	33 MPa
Tensile stress at 50% strain	3–6 MPa
Tensile strain at break	700%


Moisture Vapor Transmission Rate (MVTR) and Total Fluid Handling Capacity (TFHC)

(Method: MVTR = DIN EN 13726-2 Part 3.2

TFHC = DIN EN 13726-1 Part 3.3)



Product	Thickness	MVTR g/(day · m ²)	TFHC g/(day · m ²)
Baymedix® A on VPT 9101 T film	55 µm 25 µm	~1,800	~3,400
Baymedix® A on VPT 9101 T film	~100 µm 25 µm	~1,600	~2,800



Product	Thickness	MVTR g/(day · m ²)	TFHC g/(day · m ²)
VPT 9101 T film	25 µm	~1900	~3900



Product	Thickness	MVTR g/(day · m ²)	TFHC g/(day · m ²)
Baymedix® A	55 µm	~3,200	~10,000
Baymedix® A	~100 µm	~2,500	~6,000

The combination of the hydrophilic Baymedix® A and Platilon® film shows excellent breathability.



Covestro Deutschland AG
Business Unit Coatings,
Adhesives and Specialties
Medical
D-51365 Leverkusen
Germany

www.baymedix.com
medical@covestro.com

The manner in which you use and the purpose to which you put and utilize our products, technical assistance and information (whether verbal, written or by way of production evaluations), including any suggested formulations and recommendations, are beyond our control. Therefore, it is imperative that you test our products, technical assistance, information and recommendations to determine to your own satisfaction whether our products, technical assistance and information are suitable for your intended uses and applications. This application-specific analysis must at least include testing to determine suitability from a technical as well as health, safety, and environmental standpoint. Such testing has not necessarily been done by Covestro. All information and technical assistance is given without warranty or guarantee and is subject to change without notice. Any statement or recommendation not contained herein is unauthorized and shall not bind us. Nothing herein shall be construed as a recommendation to use any product in conflict with any claim of any patent relative to any material or its use. No license is implied or in fact granted under the claims of any patent. For more information on Covestro products in medical applications, please request from your sales support contact our guidance document: Guidance on Use of Covestro Products in a Medical Application.

Edition: 2017 · Order-No.: COV00083543 · Printed in Germany